

Meisinger

since
1888

 **MADE IN GERMANY**



Distributed by
Geistlich

MEI-BMPBA Master-Pin-Control Basic



Black Edition

The Master-Pin-Control system was developed specifically for the fixation of resorbable, non-resorbable and titanium membranes. Thanks to the sharp tip and the stable design, the pins can be placed precisely in cortical bone. The pins have an additional mini-thread and are a hybrid of screw and pin. The thread allows easy and safe removal after the successful healing period.

Product Code
MEI-BMPBA

Fixation of Resorbable Membranes, Non-Resorbable Membranes, and Titanium Mesh

Complete Tack Kit Includes: 10x Pins, 1x Autoclavable Organiser, 2x 1.8mm Initial Burs, 2x 0.8mm Twist Drills, 2x 0.6mm Twist Drills, 1x Screwdriver for Pin, 1x Pin holder, 1x Fixing holder

MEI-BMP00 Master-Pin-Control



Black Edition

Master-Pin-Control includes the same extensive set of instruments with for safe fixation of membranes as the basic kit (BMPBA) and comes with 25 additional pins.

Product Code
MEI-BMP00

Fixation of Resorbable Membranes, Non-Resorbable Membranes, and Titanium Mesh


Complete Tack Kit Includes: 35x Pins, 1x Autoclavable Organiser, 2x 1.8mm Initial Burs, 2x 0.8mm Twist Drills, 2x 0.6mm Twist Drills, 1x Screwdriver for Pin, 1x Pin holder, 1x Fixing holder




Occlusal view of severely atrophied posterior mandibular ridge




The included twist drills allow for cortical perforation to enhance and promote bleeding



Application of Geistlich Bio-Oss® and autogenous bone chips



Using the pin holder, the hybrid pins are applied into the cortical bone



The pins allow for easy insertion and safe fixation of Geistlich Bio-Gide® membrane, thus immobilising the bone graft material

Screwdriver for Pin




Fig. MEI-MP11

illustrated 1:1 (Length mm: 35.0)

Pin holder



Fig. MEI-MP12

illustrated 1:2 (Length mm: 170.0)

Fixing holder



Fig. MEI-MP14

illustrated 1:2 (Length mm: 165.0)

Initial burs
5 PK

Twist drills
2 PK

Pins
10 PK




Fig.	MEI-186RF	MEI-203RF	MEI-203RF	MEI-MP10*
Shank ¹	204	205	205	-
Size ²	018	006	008	-
Length mm	12.0	7.0	7.0	3.65
→Δ←	-	-	-	0.87
→└┘←	1.8	0.6	0.8	0.95
Head-Ø mm	-	-	-	2.51

¹ 204=RA, 205=RA L

² Largest working part diameter in 1/10

→Δ← Minimal diameter

→└┘← External diameter

Geistlich Pharma Australia & NZ

Phone AU: 1800 776 326

Phone NZ: 0800 500 043

Fax AU: 1800 709 698

Fax NZ: 0800 500 044

Email AU: info@geistlich.com.au

Email NZ: info@geistlich.co.nz