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- 1) Trullenque-Eriksson, Anna, et al. "Furcation involvement and tooth loss: A registry-based retrospective cohort study." Journal of Clinical Periodontology 50.3 (2023) : 339-347.
- 2) Jepsen, Søren, et al. "Regenerative surgical treatment of furcation defects: A systematic review and Bayesian network meta-analysis of randomized clinical trials." Journal of Clinical Periodontology 47 (2020) : 352-374.
- 3) ザ・クインテッセンス2024年5月号
- 4) Hamp, Sven-Erik, Sture Nyman, and Jan Lindhe. "Periodontal treatment of multi-rooted teeth. Results after 5 years." Journal of clinical periodontology 2.3 (1975) : 126-135.
- 5) Reynolds, Mark A., et al. "Periodontal regeneration—intrabony defects: practical applications from the AAP regeneration workshop." Clinical advances in periodontics 5.1 (2015) : 21-29.
- 6) Sanz, Mariano, et al. "Treatment of stage I-III periodontitis—The EFP S3 level clinical practice guideline." Journal of clinical periodontology 47 (2020) : 4-60.
- 7) Huang, Ren-Yeong, et al. "Mandibular disto-lingual root: a consideration in periodontal therapy." Journal of periodontology 78.8 (2007) : 1485-1490.
- 8) Casarin, Renato Corrêa Viana, et al. "Enamel matrix derivative proteins for the treatment of proximal class II furcation involvements: A prospective 24-month randomized clinical trial." Journal of clinical periodontology 37.12 (2010) : 1100-1109.
- 9) Walter, Clemens, et al. "Cone beam computed tomography (CBCT) for diagnosis and treatment planning in periodontology: systematic review update." Clinical oral investigations 24 (2020) : 2943-2958.
- 10) Sculean, Anton, et al. "Biomaterials for promoting periodontal regeneration in human intrabony defects: a systematic review." Periodontology 2000 68.1 (2015) : 182-216.
- 11) Jakse, Norbert, et al. "Primary wound healing after lower third molar surgery: evaluation of 2 different flap designs." Oral Surgery, Oral Medicine, Oral Pathology, Oral Radiology, and Endodontology 93.1 (2002) : 7-12.
- 12) Karapataki, Sofia, et al. "Healing following GTR treatment of intrabony defects distal to mandibular 2nd molars using resorbable and non-resorbable barriers." Journal of clinical periodontology 27.5 (2000) : 333-340.
- 13) Horta, Ricardo, et al. "The retromolar trigone: anatomy, cancer treatment modalities, reconstruction, and a classification system." Journal of Craniofacial Surgery 27.4 (2016) : 1070-1076.
- 14) Han, Ziyao, et al. "The keratinized mucosa width at partially edentulous molar sites and its associated factors: a pilot study." BMC Oral Health 22.1 (2022) : 630.
- 15) Kawanabe, Dai, and Ryutaro Kuraji. "The Last Molar-Entire Pad Preservation Technique (L-EPPT) as a Regenerative Surgical Approach for Combined Distal Intrabony Defects and Furcation Involvement." International Journal of Periodontics & Restorative Dentistry 45.2 (2025).
- 16) ザ・クインテッセンス2025年9月号
- 17) Yoichi T. et al. Development of Er:YAG Laser-Assisted Bone Regenerative Therapy (Er-LBRT) in Periodontal and Peri-Implant Therapy. JJSLSM Vol.40 No.1 (2019)

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