

TRI®
Performance Line



Narrow, Bone- and Tissue-Level

The TRI® Performance Line is incorporated in the Bone-Level (TRI®-Narrow, TRI®-Vent) and Tissue-Level (TRI®-Octa) implant lines and ensures optimal esthetic results both in the anterior and posterior area.



Lean & intelligent implant concept.

The TRI® Performance Line combines the pure simplicity of one connection with the esthetic flexibility of the TRI® Soft Tissue Concept. A unique lean, and yet intelligent blend for success in daily practice.



Surgical Kit

The TRI® Performance Line is characterized by one unique surgical kit for bone- and tissue-level line. Its intuitive and elegant design prioritizes efficiency and ease of use.

Implant System

TRI®-Vent Implants - Includes Surgical Cover Screw TVSCS



Implant Ø	Platform	Connection	6.5mmL	8mmL	10mmL	11.5mmL	13mmL	16mmL
3.75 mm	3.5 mmD	2.5 mmD Internal Hexagon		TV37M08	TV37M10	TV37M11	TV37M13	TV37M16
4.1 mm			TV41M06	TV41M08	TV41M10	TV41M11	TV41M13	TV41M16
4.7 mm			TV47M06	TV47M08	TV47M10	TV47M11	TV47M13	TV47M16

TRI®-Narrow Implants - Includes Surgical Cover Screw TNSCS



Implant Ø	Platform	Connection			10mmL	11mmL	13mmL	16mmL
3.3 mm	3.0 mmD	2.1 mmD Internal Hexagon			TN33M10	TN33M11	TN33M13	TN33M16

TRI®-Octa Implants - Includes Surgical Cover Screw TOSCS



Implant Ø	Platform	Connection	6.5mmL	8mmL	10mmL	11.5mmL	13mmL	
3.75 mm	4.8 mmD 1.8mm Shoulder	3.1 mmD Internal Octagon		TO37M08	TO37M10	TO37M11	TO37M13	
4.1 mm			TO41M06	TO41M08	TO41M10	TO41M11	TO41M13	
4.7 mm			TO47M06	TO47M08	TO47M10	TO47M11	TO47M13	

TRI®-Vent Prosthetics
3.5 mmD Implant Platform

HEALING COLLARS

Tissue healing - Healing Collars, sterile


TVHC35-30


TVHC50-60

Product Code		Ø	GH	Material	Qty
TVHC35-30	Healing Collar	3.5 mm	3.0 mm	Ti-6Al-4V	1
TVHC35-45	Healing Collar	3.5 mm	4.5 mm	Ti-6Al-4V	1
TVHC35-60	Healing Collar	3.5 mm	6.0 mm	Ti-6Al-4V	1
TVHC50-30	Healing Collar	5.0 mm	3.0 mm	Ti-6Al-4V	1
TVHC50-45	Healing Collar	5.0 mm	4.5 mm	Ti-6Al-4V	1
TVHC50-60	Healing Collar	5.0 mm	6.0 mm	Ti-6Al-4V	1

Surgical cover screws - sterile


TVSCS

Product Code		Ø	GH	Material	Qty
TVSCS	Surgical Cover Screw	3.5 mm		Ti-6Al-4V	1

IMPRESSION TRANSFER COMPONENTS, IMPLANT ANALOGS

Direct transfer for open-tray impressions TRI® 3D-Touch Impression - includes screw DTRS-TV05


TV05-50


DTRS-TV05

Product Code		Ø	GH	Material	Qty
TV05-35	Direct Transfer Component	3.5 mm		Ti-6Al-4V	1
TV05-50	Direct Transfer Component	5.0 mm		Ti-6Al-4V	1
DTRS-TV05	Direct Transfer Retaining Screw			Ti-6Al-4V	1

Indirect transfer/straight abutment for closed-tray impressions
TRI® 2in1 Impression Abutments - includes screws RS-TV10 and RS-TV10-ITC


TV10-35-ITC


TV10-50-ITC

Product Code		Ø	GH	Material	Qty
TV10-35-ITC	Indirect Transfer Component	3.5 mm		Ti-6Al-4V	1
TV10-50-ITC	Indirect Transfer Component	5.0 mm		Ti-6Al-4V	1
RS-TV10-ITC	Indirect Transfer Component Impression Screw			POM	5
RS-TV10	Replacement Retaining Screw			Ti-6Al-4V	1
RS-TV10-Lab	Replacement Retaining Screw - Lab Use - Green			Ti-6Al-4V	1

GH = Gingiva Height

Implant Analog


TVA

Product Code		Material	Qty
TVA-1	Digital Implant Analog	Steel	1
TVA-5	Digital Implant Analog (5 Pack)	Steel	5

COMPONENTS FOR CEMENT-RETAINED RESTORATIONS

Plastic Temporary Abutments - includes screw RS-TV10/RS-TV20


TV10-PTA


TV10-EHA-70


TV10-EHA-50

Product Code		Ø		Material	Qty
TV10-PTA	Temporary Abutment, Straight	4.5 mm		PEEK	1
RS-TV20	Replacement Retaining Screw			Ti-6Al-4V	1
RS-TV20-Lab	Replacement Retaining Screw - Lab Use - Red			Ti-6Al-4V	1
TV10-EHA-50*	Esthetic Healing Abutment	5.0 mm		PEEK	1
TV10-EHA-70*	Esthetic Healing Abutment	7.0 mm		PEEK	1
RS-TV10	Replacement Retaining Screw			Ti-6Al-4V	1
RS-TV10-Lab	Replacement Retaining Screw - Lab Use - Green			Ti-6Al-4V	1

Straight Abutments with TRI®-Friction - includes screw RS-TV10


TV10-35-F


TV10-50-F


TV10-41-F

Product Code		Ø	GH	Material	Qty
TV10-35-F	Straight Abutment with TRI®-Friction	3.5 mm		Ti-6Al-4V	1
TV10-50-F	Straight Abutment with TRI®-Friction	5.0 mm		Ti-6Al-4V	1
TV10-41-F	Contour Abutment, Straight with TRI®-Friction	4.5 mm	1.0 mm	Ti-6Al-4V	1
TV10-43-F	Contour Abutment, Straight with TRI®-Friction	4.5 mm	3.0 mm	Ti-6Al-4V	1
RS-TV10	Replacement Retaining Screw			Ti-6Al-4V	1
RS-TV10-Lab	Replacement Retaining Screw - Lab Use - Green			Ti-6Al-4V	1

Angled Abutments with TRI®-Friction - includes screws RS-TV10/RS-TV20


TV20-43-F


TV20-50-F

Product Code		Ø	GH	Material	Qty
TV20-50-F	Angled Abutment, Standard - 20° with TRI®-Friction	5.0 mm		Ti-6Al-4V	1
RS-TV20	Replacement Retaining Screw			Ti-6Al-4V	1
RS-TV20-Lab	Replacement Retaining Screw - Lab Use - Red			Ti-6Al-4V	1
TV20-41-F	Contour Abutment, Angled - 17° with TRI®-Friction	4.5 mm	1.0 mm	Ti-6Al-4V	1
TV20-43-F	Contour Abutment, Angled - 17° with TRI®-Friction	4.5 mm	3.0 mm	Ti-6Al-4V	1
RS-TV10	Replacement Retaining Screw			Ti-6Al-4V	1
RS-TV10-Lab	Replacement Retaining Screw - Lab Use - Green			Ti-6Al-4V	1

*sterile

COMPONENTS FOR SCREW-RETAINED RESTORATIONS

Screw-retained Multi-Unit Abutments



TV40-02



TV50-30



TV50-17



TTA



MUA



PCC



DTC-SRA



DTRS-SRA

Product Code		Ø	GH	Material	Qty
TV40-01*	Screw-Retained Abutment, Straight	4.5 mm	1.0 mm	Ti-6Al-4V	1
TV40-02*	Screw-Retained Abutment, Straight	4.5 mm	2.0 mm	Ti-6Al-4V	1
TV40-04*	Screw-Retained Abutment, Straight	4.5 mm	4.0 mm	Ti-6Al-4V	1
TV40-06*	Screw-Retained Abutment, Straight	4.5 mm	6.0 mm	Ti-6Al-4V	1
TV50-17*	Screw Retained Abutment, 17° Angled	4.5 mm	2.0 mm	Ti-6Al-4V	1
TV50-30*	Screw Retained Abutment, 30° Angled	4.5 mm	3.0 mm	Ti-6Al-4V	1
HC-SRAS*	Healing Cap for Screw-Retained Abutment			Ti-6Al-4V	1
TV40-SCAN	Multi-Unit Scanbody			PEEK	1
TTA	Titanium Temporary Abutment incl. Screw RS-TTA			Ti-6Al-4V	1
RS-TTA	Replacement Retaining Screw for TTA			Ti-6Al-4V	1
RS-TTAL	Replacement Retaining Screw Long for TTA			Ti-6Al-4V	1
DTC-SRA	Direct Transfer Component for Screw-Retained Abutment incl. DTRS-SRA			Ti-6Al-4V	1
DTRS-SRA	Direct Transfer Retaining Screw			Ti-6Al-4V	1
MUA	Multi-Unit Digital Implant Analog			Ti-6Al-4V	1
MUA-5	Multi-Unit Digital Implant Analog			Ti-6Al-4V	5
RS-MUA-D-4	Replacement Retaining Screw for MUA			Ti-6Al-4V	4
PCC	Plastic Castable Coping incl. Screw RS-PCC			POM	1
RS-PCC	Replacement Retaining Screw for PCC			Ti-6Al-4V	1
RS-PCC-Lab	Replacement Retaining Screw - Lab Use - Blue			Ti-6Al-4V	1

*sterile

OVERDENTURE RESTORATIONS

TRI® Raptor Abutments - includes Caps Assortment Kit (Prod. Code 192ECE)





192ECE



144AE



044CAIN

Product Code		Ø	GH	Material	Qty
130CV.01	TRI® Raptor Abutment		1.0 mm	TiN	1
130CV.02	TRI® Raptor Abutment		2.0 mm	TiN	1
130CV.03	TRI® Raptor Abutment		3.0 mm	TiN	1
130CV.04	TRI® Raptor Abutment		4.0 mm	TiN	1
130CV.05	TRI® Raptor Abutment		5.0 mm	TiN	1
130CV.06	TRI® Raptor Abutment		6.0 mm	TiN	1
140CEV	Retentive Caps/violet (strong retention 2700g)			Nylon	4
140CET	Retentive Caps/clear (standard retention 1800g)			Nylon	4
140CER	Retentive Caps/pink (soft retention 1200g)			Nylon	4
140CEG	Retentive Caps/yellow (extra-soft retention 600g)			Nylon	4
192ECE	Caps Assortment Kit			Ti, Nylon	1
335SBC	Smart Box Kit			Ti, Nylon	1
144AE	Analog for Lab			Stainless	2
044CAIN	Transfer for Impression			Stainless	2

COMPONENTS FOR TRI+ DIGITAL CAD/CAM

CAD/CAM Abutments - includes screw RS-TV10



TV70-SCAN



TV70-20-F



TV70-07-NE



TV70-Blank-F

Product Code		Ø	GH	Material	Qty
TV70-SCAN	Scanbody for Laboratory and Intraoral Scanner			PEEK	1
TV70-07-F	Engaging Titanium Bonding-Base with TRI®-Friction	4.5 mm	0.7 mm	Ti-6Al-4V	1
TV70-20-F	Engaging Titanium Bonding-Base with TRI®-Friction	4.5 mm	2.0 mm	Ti-6Al-4V	1
TV70-07-NE	Non-Engaging Titatium Bonding-Base	4.5 mm	0.7 mm	Ti-6Al-4V	1
TV70-20-NE	Non-Engaging Titatium Bonding-Base	4.5 mm	2.0 mm	Ti-6Al-4V	1
TV70-Blank-F	Titanium Milling Blank with TRI®-Friction			Ti-6Al-4V	1
RS-TV10	Replacement Retaining Screw			Ti-6Al-4V	1
RS-TV10-Lab	Replacement Retaining Screw - Lab Use - Green			Ti-6Al-4V	1

GH = Gingiva Height

TRI®-Narrow Prosthetics
3.0 mmD Implant Platform

HEALING COLLARS

Tissue healing - Healing Collars, sterile


TNHC35-30

Product Code		Ø	GH	Material	Qty
TNHC35-30	Healing Collar	3.5 mm	3.0 mm	Ti-6Al-4V	1
TNHC35-45	Healing Collar	3.5 mm	4.5 mm	Ti-6Al-4V	1
TNHC35-60	Healing Collar	3.5 mm	6.0 mm	Ti-6Al-4V	1
TNHC45-30	Healing Collar	4.5 mm	3.0 mm	Ti-6Al-4V	1
TNHC45-50	Healing Collar	4.5 mm	5.0 mm	Ti-6Al-4V	1

Surgical cover screws, sterile


TNSCS

Product Code		Ø	GH	Material	Qty
TNSCS	Surgical Cover Screw	3.2 mm		Ti-6Al-4V	1

IMPRESSION TRANSFER COMPONENTS, IMPLANT ANALOGS

Direct transfer for open-tray impressions - TRI® 3D-Touch Impression - includes screw DTRS-TN05


TN05-35

DTRS-TN05

Product Code		Ø	GH	Material	Qty
TN05-35	Direct Transfer Component	3.5 mm		Ti-6Al-4V	1
DTRS-TN05	Direct Transfer Retaining Screw			Ti-6Al-4V	1

Indirect transfer/straight abutment for closed-tray impressions
TRI® 2in1 Impression Abutments - includes screw RS-TN10


TN 10-35-ITC

Product Code		Ø	GH	Material	Qty
TN10-35-ITC	Indirect Transfer Component	3.5 mm		Ti-6Al-4V	1
RS-TN10	Replacement Retaining Screw			Ti-6Al-4V	1

Implant Analog


TNA

Product Code		Material	Qty
TNA-1	Digital Implant Analog	Ti-6Al-4V	1
TNA-5	Digital Implant Analog (5 Pack)	Ti-6Al-4V	5

GH = Gingiva Height

COMPONENTS FOR CEMENT-RETAINED RESTORATIONS

Plastic Temporary Abutments - includes screw RS-TN20


TN10-PTA

Product Code		Ø		Material	Qty
TN10-PTA	Temporary Abutment, Straight	3.5 mm		PEEK	1
RS-TN20	Replacement Retaining Screw			Ti-6Al-4V	1
RS-TN20-Lab	Replacement Retaining Screw - Lab Use - Red			Ti-6Al-4V	1

Straight Abutments with TRI®-Friction - includes screw RS-TN10


TN10-41-F


TN10-35-F

Product Code		Ø	GH	Material	Qty
TN10-35-F	Straight Abutment with TRI®-Friction	3.5 mm		Ti-6Al-4V	1
TN10-41-F	Contour Abutment, Straight with TRI®-Friction	4.0 mm	1.0 mm	Ti-6Al-4V	1
TN10-43-F	Contour Abutment, Straight with TRI®-Friction	4.0 mm	3.0 mm	Ti-6Al-4V	1
RS-TN10	Replacement Retaining Screw			Ti-6Al-4V	1
RS-TN10-Lab	Replacement Retaining Screw - Lab Use - Green			Ti-6Al-4V	1

Angled Abutments with TRI®-Friction - includes screws RS-TN10


TN20-41-F


TN20-45-F

Product Code		Ø	GH	Material	Qty
TN20-45-F	Angled Abutment, Standard - 20° with TRI®-Friction	4.5 mm		Ti-6Al-4V	1
TN20-41-F	Contour Abutment, Angled - 17° with TRI®-Friction	4.0 mm	1.0 mm	Ti-6Al-4V	1
TN20-43-F	Contour Abutment, Angled - 17° with TRI®-Friction	4.0 mm	3.0 mm	Ti-6Al-4V	1
RS-TN10	Replacement Retaining Screw			Ti-6Al-4V	1
RS-TN10-Lab	Replacement Retaining Screw - Lab Use - Green			Ti-6Al-4V	1

COMPONENTS FOR TRI+ DIGITAL CAD/CAM

CAD/CAM Abutments - includes screw RS-TN10


TN70-SCAN


TN70-07-F


TN70-07-NE


TN70-Blank-F

Product Code		Ø	GH	Material	Qty
TN70-SCAN	Scanbody for Laboratory and Intraoral Scanner			PEEK	1
TN70-07-F	Engaging Titanium Bonding-Base with TRI®-Friction	4.2 mm	0.7 mm	Ti-6Al-4V	1
TN70-07-NE	Non-Engaging Titatium Bonding-Base	4.2 mm	0.7 mm	Ti-6Al-4V	1
TN70-Blank-F	Titanium Milling Blank with TRI®-Friction			Ti-6Al-4V	1
RS-TN10	Replacement Retaining Screw			Ti-6Al-4V	1
RS-TN10-Lab	Replacement Retaining Screw - Lab Use - Green			Ti-6Al-4V	1

COMPONENTS FOR SCREW-RETAINED RESTORATIONS

Screw-retained Multi-Unit Abutments


TN40-02


TTA


MUA


PCC


DTC-SRA


DTRS-SRA

Product Code		Ø	GH	Material	Qty
TN40-02*	Screw-Retained Abutment, Straight	4.5 mm	2.0 mm	Ti-6Al-4V	1
TN40-04*	Screw-Retained Abutment, Straight	4.5 mm	4.0 mm	Ti-6Al-4V	1
HC-SRAS*	Healing Cap for Screw-Retained Abutment			Ti-6Al-4V	1
TV40-SCAN	Multi-Unit Scanbody			PEEK	1
TTA	Titanium Abutment incl. Screw RS-TTA			Ti-6Al-4V	1
RS-TTA	Replacement Retaining Screw for TTA			Ti-6Al-4V	1
RS-TTAL	Replacement Retaining Screw Long for TTA			Ti-6Al-4V	1
DTC-SRA	Direct Transfer Component for Screw-Retained Abutment incl. DTRS-SRA			Ti-6Al-4V	1
DTRS-SRA	Direct Transfer Retaining Screw			Ti-6Al-4V	1
MUA	Multi-Unit Digital Implant Analog			Ti-6Al-4V	1
MUA-5	Multi-Unit Digital Implant Analog			Ti-6Al-4V	5
RS-MUA-D-4	Replacement Retaining Screw for MUA			Ti-6Al-4V	4
PCC	Plastic Castable Coping incl. Screw RS-PCC			POM	1
RS-PCC	Replacement Retaining Screw for PCC			Ti-6Al-4V	1
RS-PCC-Lab	Replacement Retaining Screw - Lab Use - Blue			Ti-6Al-4V	1

OVERDENTURE RESTORATIONS

TRI® Raptor Abutments - includes Caps Assortment Kit (Prod. Code 192ECE)




192ECE


144AE


044CAIN

Product Code		Ø	GH	Material	Qty
130TN.01	TRI® Raptor Abutment		1.0 mm	TiN	1
130TN.02	TRI® Raptor Abutment		2.0 mm	TiN	1
130TN.03	TRI® Raptor Abutment		3.0 mm	TiN	1
130TN.04	TRI® Raptor Abutment		4.0 mm	TiN	1
130TN.05	TRI® Raptor Abutment		5.0 mm	TiN	1
140CEV	Retentive Caps/violet (strong retention 2700g)			Nylon	4
140CET	Retentive Caps/clear (standard retention 1800g)			Nylon	4
140CER	Retentive Caps/pink (soft retention 1200g)			Nylon	4
140CEG	Retentive Caps/yellow (extra-soft retention 600g)			Nylon	4
192ECE	Caps Assortment Kit			Ti, Nylon	1
335SBC	Smart Box Kit			Ti, Nylon	1
144AE	Analog for Lab			Stainless	2
044CAIN	Transfer for Impression			Stainless	2

*sterile

TRI®-Octa Prosthetics
4.8 mmD Implant Platform

HEALING COLLARS

Tissue healing - Healing Collars, sterile


TOHC51-15


TOHC51-30

Product Code		Ø	GH	Material	Qty
TOSCS	Surgical Cover Screw	3.5 mm		Ti-6Al-4V	1
TOHC51-15	Healing Collar	5.1 mm	1.5 mm	Ti-6Al-4V	1
TOHC51-30	Healing Collar	5.1 mm	3.0 mm	Ti-6Al-4V	1

IMPRESSION TRANSFER COMPONENTS, IMPLANT ANALOGS

Direct transfer for open-tray impressions - TRI® 3D-Touch Impression - includes screw DTRS-TO05


TO05-51


DTRS-TO05

Product Code		Ø	GH	Material	Qty
TO05-51	Direct Transfer Component	5.1 mm		Ti-6Al-4V	1
DTRS-TO05	Direct Transfer Retaining Screw			Ti-6Al-4V	1

Indirect transfer/straight abutment for closed-tray impressions

TRI® 2in1 Impression Abutments - includes screw RS-TO


TO10-51-ITC

Product Code		Ø	GH	Material	Qty
TO10-51-ITC	Indirect Transfer Component	5.1 mm		Ti-6Al-4V	1
RS-TO	Replacement Retaining Screw			Ti-6Al-4V	1

Implant Analog


TOA

Product Code		Material	Qty
TOA-1	Digital Implant Analog	Ti-6Al-4V	1
TOA-5	Digital Implant Analog (5 Pack)	Ti-6Al-4V	5

GH = Gingiva Height

COMPONENTS FOR CEMENT-RETAINED RESTORATIONS

Plastic Temporary Abutments - includes screw RS-TO



TO10-PTA

Product Code		Ø		Material	Qty
TO10-PTA	Temporary Abutment, Straight	5.1 mm		PEEK	1
RS-TO	Replacement Retaining Screw			Ti-6Al-4V	1

Straight Abutments with TRI® Friction - includes screw RS-TO



TO10-35-F



TO10-51-F

Product Code		Ø	GH	Material	Qty
TO10-35-F	Straight Abutment with TRI®-Friction	3.5 mm		Ti-6Al-4V	1
TO10-51-F	Straight Abutment with TRI®-Friction	5.1 mm		Ti-6Al-4V	1
RS-TO	Replacement Retaining Screw			Ti-6Al-4V	1

Angled Abutments with TRI®-Friction - includes screws RS-TO



TO20-51-F

Product Code		Ø	GH	Material	Qty
TO20-51-F	Angled Abutment, Standard - 20° with TRI®-Friction	5.1 mm		Ti-6Al-4V	1
RS-TO	Replacement Retaining Screw			Ti-6Al-4V	1

COMPONENTS FOR SCREW-RETAINED RESTORATIONS

Screw-retained Multi-Unit Abutments



TO40-16



TO40-30



TTA



MUA



PCC



DTC-SRA



DTRS-SRA

Product Code		Ø	GH	Material	Qty
TO40-16*	Screw-Retained Abutment, Straight	4.5 mm	1.6 mm	Ti-6Al-4V	1
TO40-30*	Screw-Retained Abutment, Straight	4.5 mm	3.0 mm	Ti-6Al-4V	1
HC-SRAS*	Healing Cap for Screw-Retained Abutment			Ti-6Al-4V	1
TV40-SCAN	Multi-Unit Scanbody			PEEK	1
TTA	Titanium Abutment incl. Screw RS-TTA			Ti-6Al-4V	1
RS-TTA	Replacement Retaining Screw for TTA			Ti-6Al-4V	1
RS-TTAL	Replacement Retaining Screw Long for TTA			Ti-6Al-4V	1
DTC-SRA	Direct Transfer Component for Screw-Retained Abutment incl. DTRS-SRA			Ti-6Al-4V	1
DTRS-SRA	Direct Transfer Retaining Screw			Ti-6Al-4V	1
MUA	Multi-Unit Digital Implant Analog			Ti-6Al-4V	1
MUA-5	Multi-Unit Digital Implant Analog			Ti-6Al-4V	5
RS-MUA-D-4	Replacement Retaining Screw for MUA			Ti-6Al-4V	4
PCC	Plastic Castable Coping incl. Screw RS-PCC			POM	1
RS-PCC	Replacement Retaining Screw for PCC			Ti-6Al-4V	1
RS-PCC-Lab	Replacement Retaining Screw - Lab Use - Blue			Ti-6Al-4V	1

*sterile

OVERDENTURE RESTORATIONS

TRI® Raptor Abutments - includes Caps Assortment Kit (Prod. Code 192ECE)





192ECE



144AE



044CAIN

Product Code		Ø	GH	Material	Qty
130BFT.00	TRI® Raptor Abutment		0.0 mm	TiN	1
130BFT.01	TRI® Raptor Abutment		1.0 mm	TiN	1
130BFT.02	TRI® Raptor Abutment		2.0 mm	TiN	1
130BFT.03	TRI® Raptor Abutment		3.0 mm	TiN	1
130BFT.04	TRI® Raptor Abutment		4.0 mm	TiN	1
130BFT.05	TRI® Raptor Abutment		5.0 mm	TiN	1
140CEV	Retentive Caps/violet (strong retention 2700g)			Nylon	4
140CET	Retentive Caps/clear (standard retention 1800g)			Nylon	4
140CER	Retentive Caps/pink (soft retention 1200g)			Nylon	4
140CEG	Retentive Caps/yellow (extra-soft retention 600g)			Nylon	4
192ECE	Caps Assortment Kit			Ti, Nylon	1
335SBC	Smart Box Kit			Ti, Nylon	1
144AE	Analog for Lab			Stainless	2
044CAIN	Transfer for Impression			Stainless	2

COMPONENTS FOR TRI+ DIGITAL CAD/CAM

CAD/CAM Abutments - includes screw RS-TO



TO70-SCAN



TO70-07-F



TO70-07-NE




TO70-Blank-F


Product Code		Ø	GH	Material	Qty
TO70-SCAN	Scanbody for Laboratory and Intraoral Scanner			PEEK	1
TO70-07-F	Engaging Titanium Bonding-Base with TRI®-Friction	5.1 mm	0.7 mm	Ti-6Al-4V	1
TO70-07-NE	Non-Engaging Titatium Bonding-Base	5.1 mm	0.7 mm	Ti-6Al-4V	1
TO70-Blank-F	Titanium Milling Blank with TRI®-Friction			Ti-6Al-4V	1
RS-TO	Replacement Retaining Screw			Ti-6Al-4V	1

GH = Gingiva Height

Drill & Implant Lengths



Torque Requirements




Drill Sequence


3.3 mmD



3.75 mmD



4.1 mmD



4.7 mmD

